



PRECISION
FOR A
HEALTHIER
TOMORROW

ORTHODONTIC INSTRUMENTS

Precision Tools for Confident Smiles



PRECISION
Engineered
for Excellence



QUALITY
Trusted by
Professionals



BETTER SMILES
Supporting
Healthier Futures



Qbro Industries



Phone/WhatsApp
+4917620793624
+923007705082
+923039359066



info@qbroindustries.com

DENTAL SOLUTIONS
FOR A BRIGHTER SMILE

ORTHODONTIC INSTRUMENTS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان



QB-0400
Height Gauge 0.18
Aluminum



QB-0401
Height Gauge 0.22
Aluminum



QB-0402
Height Gauge 0.18
Stainless Steel



QB-0403
Height Gauge 0.22
Stainless Steel



QB-0404
Moveable
Height Gauge



QB-0405
Endo Finger
Ruler Span
Measure Scale



QB-0406
Endo Finger
Ruler Span
Measure Scale



QB-0407
Endo Finger
Ruler Span
Measure Scale



QB-0408
Positive
Postioner



QB-0409
Strip Saw



QB-0410
Bone Bracket Gauge



QB-0411
Endo Finger Ruler
Span Measure Scale



QB-0412
Endo Finger Ruler
Span Measure Scale

ORTHODONTIC INSTRUMENTS

- Щипцы для ортодонтии и протезирования
- Pincettes pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان



QB-0413
Measuring Loop
angled, length
175 mm,
measuring range
0 - 40 mm



QB-0414
Measuring Loop
straight length
175 mm,
measuring range
0 - 40 mm



QB-0415
Measuring loop, Castroviejo
angled, length 175 mm,
measuring range 0 - 40 mm,
with window



QB-0416
Measuring Loop
straight length 90 mm,
measuring range
0 - 40 mm,

QB-0417
Castroviejo Calipers
20mm
90mm



QB-0418
Measuring Loop
angled, length 90 mm,
measuring range
0 - 40 mm,

QB-0419
Castroviejo Calipers
20mm, 45°
90mm



QB-0420
Castroviejo Calipers
20mm
90mm



QB-0421
Fig. 1
Iwanson
spring caliper
For metal and
porcelain.



QB-0422
Fig. 2
Iwanson
spring caliper
For wax.

ORTHODONTIC INSTRUMENTS

- Щипцы для ортодонтии и протезирования
- Pincettes pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان



QB-0423
Boley Gauge



QB-0424
Vernier Gauge



QB-0425
Aluminium Calliper
up to 200mm



QB-0426
Boley Gauge



QB-0427
Delux Boley
Gauge



QB-0428
Beerendonk
130mm
Mini Gauge



QB-0429
Compasso
De Willis



QB-0430
Dentaum
Ruler



QB-0431
Ruler Wedge



QB-0432
Ruler

ORTHODONTIC INSTRUMENTS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان

Bracket Holder Tweezers



QB-0433
Deluxe Bracket
Holder Tweezers



QB-0434
Deluxe Bracket
Holder Tweezers



QB-0435
Direct Bond Bracket
Holder Tweezers



QB-0436
Bracket Tweezers
Slot Aligner



QB-0437
Direct Bond
Bracket Holder
Tweezers, Small Tip



QB-0438
Exact Bracket
Holder Tweezers



QB-0439
Bracket Placing
Tweezers



QB-0440
Easy Access
Tweezers



QB-0441
Bracket Placer
Tweezers
(Extra Fine Tip)



QB-0442
Deluxe Bracket
Holder Tweezers
135mm



QB-0443
Easy Eccess
Tweezers



QB-0444
Slipper
Tweezers



QB-0445
Direct Bonding
Bracket Holders
135mm



QB-0446
Direct Bonding
Bracket Holders
136mm



QB-0447
Direct Bonding
Bracket Holders
135mm



QB-0448
Deluxe Bracket
Holder



QB-0449
Direct Bonding
Bracket Holders
140mm

ORTHODONTIC INSTRUMENTS

- Щипці для ортодонції та протезування
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان



QB-0450
Band seater
auto cleavable



QB-0451
Band seater
auto cleavable
Triangle Shape



QB-0452
Band seater
auto cleavable
Square Shape



QB-0453
Band seater
auto cleavable
Round Shape



QB-0454
Mershon Band
Pusher
140mm



QB-0455
Band seating &
removing
instruments



QB-0456
Ortho Light double
ended
ligature director
scaler



QB-0457
Ortho Light Band
Pusher Scaler



QB-0458
Band Seating &
Removing
Silicon Grip



QB-0459
Ligature Instrument
Silicon Grip



QB-0460
Elastic Instrument
Silicon Grip



QB-0461
Distal Bender
Silicon Grip



QB-0462
Distal Bender
Silicon Grip



QB-0463
Distal Bender
Silicon Grip

ORTHODONTIC INSTRUMENTS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان



QB-0464
Ortho light scaler double ended



QB-0465
Ortho light scaler double ended



QB-0466
Ortho light scaler double ended



QB-0467
Twiel-on elastic module instrument



QB-0468
Twiel-on elastic module instrument



QB-0469
Schure Band Seater



QB-0470
Ortho light scaler



QB-0471
Ortho light scaler



QB-0472
Ortho light scaler heavy



QB-0473
Ortho instrument



QB-0474
Ortho instrument

ORTHODONTIC INSTRUMENTS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان



QB-0475
Ortho Light double ended
ligature director



QB-0476
Ortho Light double ended ligature scaler director tivate coil spring
are remove modules with one instruments



QB-0477
Ortho band pusher scaler



QB-0478
Ortho ligature director single ended



QB-0479
Ortho Light double ended ligature
director scaler



QB-0480
Ortho Light double ended ligature scaler/director
tivate coil spring are remove modules with one instruments



QB-0481
Ortho Light double ended ligature director



QB-0482
Ortholite Ligature Director Scaler



QB-0483
Ortho light ligature director double ended



QB-0484
Double ended cement spatula



QB-0485
Ortho light scaler single ended

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان

BOX JOINT PLIERS

QB-0486
Universal Dolphin Style
Wire Bending & cutting pliers



QB-0487 Band Removing Plier
The curved blade of the band removing plier can strip easily and gently remove directly to the tooth. Length 130mm



QB-0488
Tweed Loop forming Pliers



QB-0489
How Plier curved
Universally usable tool for the multiband technology for wires up to .028". Length 130mm



QB-0490 How Plier
Universally usable tool for the multiband technology for wires up to .028". Length 135mm



QB-0491
Nance Loop forming Pliers



QB-0492
Bracket Removing Plier curved
Bracket Removing Plier curved for easy removal of brackets in the posterior region. Length 130mm



QB-0493
Young Wire Bending & Loop Forming Pliers



QB-0494
Wire Bending & Loop Forming Pliers



QB-0495
Adams, small wire bending and loop forming Pliers



QB-0496
Bracket Removing Plier Round Box
Bracket removing plier for comfortable and quick removal of the brackets. Length 135mm



QB-0497
Weingart Plier medium
Weingart plier that are ideal due to their curved tip for setting sheets and is capable of bending the ends. Length 145mm



PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان

BOX JOINT PLIERS

QB-0498 Adams Plier Round Box

Adams pliers for bending wire arches and bends right angles for wires up to .022", .025 tip to beak. Length 125mm



QB-0499 Falcon Beak Contouring Pliers.



QB-0500 wires up to 0.020" / 0.7mm.



QB-0501 Johnson Contouring Pliers



QB-0502 Aderer Large wire bending and loop forming pliers



QB-0503 Tweed / angle rectangular arch forming & Torquing pliers



QB-0504 Adams Plier

Adams pliers for bending wire arches and bends right angles for wires up to .022", .025 tip to beak. Length 125mm



QB-0505 Tweed Loop Formin Pliers



QB-0506 Bird Beak Plier Standard Bird Beak Wire bending plier without notches, which is suitable as a universal utility tool. Length 125mm



QB-0507 Mini Ligature cutters



QB-0508 Small Fluch-cut distal and cutting Pliers with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips



QB-0509 Large Fluch-cut distal and cutting Pliers with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips



PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان

BOX JOINT PLIERS

QB-0510

Tweed / Angel

1.78mm (.070") wide tips arch forming and torquing pliers



QB-0516

Contouring Plier



QB-0511

Tweed / Angel

1.78mm (.070") wide tips arch forming and torquing pliers



QB-0517

Bands & Crowns Plier
in order to skirt



QB-0512

Lingual / Palatal arch
forming pliers



QB-0518

Ligature Forming Pliers



QB-0513

Aderer

Facebow loop forming Pliers



QB-0519

for double back and triple back bands in 0.30" / 0.76mm & 0.36" / 0.91mm lingual arch wires



QB-0514

Waldsavhs wire bending and cutting pliers



QB-0520

Bird Beak Wire Bending & cutting pliers



QB-0515 **Aderer Plier**

For the manufacturing of appliances and activation of wires for wires up to .030" fine tip and up to .036" mean tip. Length 130mm



QB-0521

Weingart Plier medium

Weingart plier that are ideal due to their curved tip for setting sheets and is capable of bending the ends. Length 145mm



PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان

L-KEY JOINT PLIERS



QB-0522

Distal end cutter with safety hold, mini head, long handles.
fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" x .025")
- ⊗ 0.38 mm to 0.53 mm (.0155" to .0215")
- ⊗ Max. 0.53 mm x 0.64 mm (.021" x .025")



QB-0526

Soft wire and ligature cutter, 15° angled fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



QB-0523

Flush-cut distal end cutter with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" to .025")



QB-0527

Soft wire and ligature cutter, 45° angled fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



QB-0524

Soft wire and ligature cutter, Mini fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



QB-0528

Hard wire cutter, straight Fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



QB-0525

Soft wire and ligature cutter, standard fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ⌋ .25 mm to .38 mm (.010" to .015")



QB-0529

Hard wire cutter, 15° angled Fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



QB-0530

Hard wire cutter, 15° angled Non-autoclavable

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان

L-KEY JOINT PLIERS



QB-0531
Band removing pliers



QB-0536
Weingart-mini
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



QB-0532
Bracket removing pliers angled



QB-0537
Weingart
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



QB-0533
Bracket removing pliers straight



QB-0538
How Straight
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



QB-0534
Bracket removing pliers with plastic tip



QB-0539
Angled
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



QB-0535
Adhesive removing pliers



QB-0540
Wire bending and loop forming pliers
Angled (Bird beak)

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان

L-KEY JOINT PLIERS



QB-0541

Multi-purpose pliers, bird beak,
De la rosa and cutter in one

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.56 mm (.018" x .022")



QB-0546

Stop pliers V-Bebd - 1mm

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.43 mm x 0.56 mm (.017" x .028")



QB-0542

Cinch Back / tie back
Distal end bending pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



QB-0547

Lingual / palatal arch forming plies

- 0.76 mm (.030") / 0.91 mm (.036")



QB-0543

Hammer head
Distal end bending pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



QB-0548

Torquing pliers set



QB-0544

Nance
Loop forming pliers

- 0.36 mm to 0.71 mm (.014" to .028")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



QB-0549

Tweed / angel
Rectangular arch forming and
torquing pliers

- 0.30 mm to 0.71 mm (.012" to .028")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



QB-0545

Detailing step pliers / Bayonet
bending pliers - 1mm

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



QB-0551

Surgical ball hook crimping
pliers, straight



QB-0550

Surgical ball hook crimping pliers,
angled



PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان

L-KEY JOINT PLIERS




QB-0552
Aderer Clamp three jaw wire bending plier for contouring and bending wire up to .30"/.76mm.



QB-0558
Falcon beak contouring plier.




QB-0553
Bird beak plier for bending round, square and rectangle wire up to .30"/.76mm.



QB-0559
For placement and removal of arch wires, pins and other auxiliaries.




QB-0554
Bird beak plier for light wire bending round, square and rectangular up to .016"/0.41mm.



QB-0560
For placement and removal of arch wires, pins and other auxiliaries in angle design



QB-0555
Double Rounded jaw bending pliers, hooks and wire up to .020"/0.51mm.



QB-0561
contouring Plier.



QB-0556
Adams pliers for hooks and bending wire up to .030"/0.76mm.



QB-0563
Loop forming in four length 3,4,5 & 6 mm wire capacity 0.22", 56mm.



QB-0557
Johnson pliers for crown band contouring



QB-0564
Arch bending plier design for convenient handling of square or rectangular wire, .020" / 0.51mm.



PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان

LAP & L-KEY JOINT PLIERS



QB-0565
In order to skirt
bands and crowns.



QB-0566
Wires up to
.020" / 0.7mm.



QB-0567
Round and concave plier
for loop forming and banding
wire up to .20" / .051mm.



QB-0568
Bending for round square
or rectangular wire up
to 0.30"/0.76mm



QB-0569
for double back and triple
back bands in 0.30" / .76mm
& 0.36" / .91mm lingual arch wires.



QB-0570
Detailing step plier double-sided
offset beak permits bayonet bends.
Creates right hand or left hand
step-up or stepdown at 1/4, 1/2, or
1mm increments.



QB-0571
"V" bend stop special design makes
1mm "V" bends to shorten archwires
or provide a positive stop. Excellent
for placing stops in nitium wires.
For round or rectangular wires up to .025"



QB-0572
Angle how plier efficient utility plier
with serrated beaks for a firm grip.
45 angled tips to access hard to
reach areas.



QB-0573
Weingardt plier serrated tips hold
wires firmly at a convenient working
angle. Tapered beaks fit easily
between brackets and hard to reach
areas. Our most popular utility plier.



QB-0574
Pin & Ligature Cutters
TC.



QB-0575
Mini Pin & Ligature
Cutter TC.



QB-0576
Standard Ligature
Cutter TC.



QB-0577
V Bend Stop Plier
TC.




QB-0578
Angled Direct Bond Bracket
Remover TC.



QB-0579
Cutter
Str.



Distal Cutter
14cm
QB-0580 2.5 mm
QB-0581 3.0 mm
QB-0582 3.0 mm



QB-0583
Arrow Clasp
Forming Pliers
12.0cm



QB-0584
Schwarz
12.5cm

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- Pinzetten voor orthodontie en prothese

QB-0585

Distal end Cutter, Safety Hold Plier
(Ideal for Cutting Wire) TC.



QB-0586

Hook Crimping Plier Angled
Straight TC.



QB-0587

3 Prong Clasp Adjustable
Plier TC.



QB-0588

Utility Plier Weingart TC.



QB-0589

Direct Bracket Remover
TC.



QB-0590

Straight Direct Bond Bracket
Remover TC.



QB-0591

Pin Holding Plier TC.



QB-0592

Hard wire cutters



QB-0593

Elastic Separating Plier
Round Handle



QB-0594

Elastic
Separating Plier



QB-0595 Str.
QB-0596 Cvd.
Smaha
Elastomeric
ligature
placing forceps
140 mm



QB-0597

Elastomeric
ligature placing
forceps with
groove
Falcon 140 mm



QB-0598

Elastomeric
ligature placing
forceps with
groove
Falcon 140 mm



PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pinces pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان

CUTTER, UTILITY PLIERS & WIRE FORMING PLIERS



QB-0599
Weingart Pliers



QB-0604
Adhesive
Removing Pliers



QB-0600
Slim Weingart
Pliers



QB-0605
Posterior Band Removing
Pliers - Long



QB-0601
How Pliers
Straight



QB-0606
Posterior Band Removing
Pliers - Short




QB-0602
How Pliers Offset



QB-06077
Bird Beak Pliers



QB-0603
Band seating Pliers



QB-06088
Omega Loop
Forming Pliers

Replacement tip and
hex key included.

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincettes pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان

CUTTERS




QB-0609
Hard
Wire Cutter 15




QB-0613
Mini Pin &
Ligature Cutter



QB-0610
Distal End Cutter
Flush Cut



QB-0614
Hard Wire Cutter
Straight



QB-0611
Universal Cut & Hold
Distal End Cutter



QB-0615
Micro-Mini Pin &
Ligature Cutter



QB-0612
Pin & Ligature
Cutter



QB-0616
Cutter

PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان



QB-0617
Three Jaw Pliers



QB-0621
Jarabak Pliers



QB-0618
Light Wore Bird Beak
Pliers with Groove



QB-0622
Hollow Chop
Contouring Pliers




QB-0619
Tweed loop
Forming Pliers



QB-0623
Lingual Arch
Forming Pliers



QB-0620
Arch Bending Pliers



QB-0624
Rectangular Arch
Bending Pliers

ORTHODONTICS PLIERS STANDS

- Щипцы для ортодонтии и протезирования
- Pincettes pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik

- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- أدوات تقويم الأسنان



QB-0625
Plier Stand
Stainless Steel.



QB-0628
Plier Stand
Stainless Steel.



QB-0626
Pliers Stand
Stainless Rod
(Aluminum Anodized)
Stainless Steel



QB-0629
Plier Stand
Stainless Steel.



QB-0627
Plier Stand
Stainless Steel.



QB-0630
Plier Stand
Plastic



PRECISION
FOR A
HEALTHIER
TOMORROW

Trusted Orthodontic Solutions

At Qbro Industries, we are committed to precision, quality and innovation. We design and manufacture reliable orthodontic instruments that support dental professionals worldwide in creating healthier smiles and brighter futures.

DENTAL SOLUTIONS
FOR A BRIGHTER SMILE



PRECISION

Engineered
for Excellence



QUALITY

Trusted by
Professionals



INNOVATION

Driving Better
Smiles



Qbro Industries



Phone/WhatsApp

+4917620793624

+923007705082

+923039359066



info@qbroindustries.com



www.qbroindustries.com



SCAN TO CONNECT

BETTER INSTRUMENTS. BRIGHTER SMILES.