



# Qbro

— INDUSTRIES —

INSTRUMENTS  
FOR A HEALTHIER  
BRIGHTER SMILE

BETTER  
TOOLS  
BRIGHTER  
SMILES

## DIAGNOSTIC, PERIODONTAL & RESTORATIVE

*Trusted  
by Professionals  
Worldwide*

### — INSTRUMENTS —

Advanced Instruments for Accurate Diagnosis  
and Reliable Restorations



**PRECISION**  
IN EVERY DETAIL



**QUALITY**  
YOU CAN TRUST



**INNOVATION**  
FOR A BRIGHTER  
TOMORROW



GLOBAL  
PARTNER  
IN DENTAL  
EXCELLENCE

#### Qbro Industries




Phone/WhatsApp:  
+4917620793624  
+923007705082  
+923039359066



Email: [info@qbroindustries.com](mailto:info@qbroindustries.com)



QUALITY CARE  
WORLDWIDE



## MOUTH MIRRORS

Plain



QB-6600  
Fig. 3  
Ø20mm

QB-6601  
Fig. 4  
Ø22mm

QB-6602  
Fig. 5  
Ø24mm

QB-6603  
Fig. 6  
Ø26mm



Cone-Socket Plain



QB-6604  
Fig. 3  
Ø20mm

QB-6605  
Fig. 4  
Ø22mm

QB-6606  
Fig. 5  
Ø24mm

Magnifying



QB-6607  
Fig. 3  
Ø20 mm

QB-6608  
Fig. 4  
Ø22 mm

QB-6609  
Fig. 5  
Ø24 mm



Cone-Socket Magnifying



QB-6610  
Fig. 3  
Ø20mm

QB-6611  
Fig. 4  
Ø22mm

QB-6612  
Fig. 5  
Ø24mm

Rhodium



QB-6613  
Fig. 3  
Ø20mm

QB-6614  
Fig. 4  
Ø22mm

QB-6615  
Fig. 5  
Ø24 mm



Cone-Socket Rhodium



QB-6616  
Fig. 3  
Ø20mm

QB-6617  
Fig. 4  
Ø22mm

QB-6618  
Fig. 5  
Ø24mm

Front Surface



QB-6619  
Fig. 3  
Ø20mm

QB-6620  
Fig. 4  
Ø22mm

QB-6621  
Fig. 5  
Ø24mm

QB-6622  
Fig. 6  
Ø26mm



Cone-Socket Front Surface



QB-6623  
Fig. 3  
Ø20mm

QB-6624  
Fig. 4  
Ø22mm

QB-6625  
Fig. 5  
Ø24mm

## HANDLES FOR MOUTH MIRRORS



QB-6626  
Hollow Mirror  
Handle Round  
Ø 10 mm



QB-6627  
Hollow Mirror  
Hexagonal  
Handle  
Ø 8 mm



QB-6628  
Hollow  
Ø 8 mm



QB-6629  
Hollow Mirror  
Hexagonal  
Handle  
Ø 10 mm



QB-6630  
Hollow Mirror  
Handle  
Ø 8 mm



QB-6631  
Hollow Mirror  
Handle  
Ø 10 mm



QB-6632  
Mirror Handle  
Solid  
Ø 6 mm



QB-6633  
Mirror Handle  
Solid  
Ø 6 mm



QB-6634  
Mirror Handle  
Solid Octagonal  
Ø 6 mm



QB-6635  
Mirror Handle  
Stem  
with Rule  
35mm



QB-6636  
Mirror Handle  
Solid  
Stem with  
Rule 35mm



QB-6637  
Cone-Socket  
Solid Octagonal  
Ø 6 mm



QB-6638  
Cone-socket  
Solid  
Ø 6 mm



QB-6639  
Cone-socket  
Hexagonal  
Handle  
Ø 6 mm

Cone-Socket Mirror Handle

## EXPLORERS

Examination kit Set of 3 Pieces



QB-6640  
Probe

QB-6641  
Mirror

QB-6642  
College Forcep  
15cm



Examination Set of 4 Pieces



QB-6643  
Exp,23/17-A

QB-6644  
Excavators

QB-6645  
Tweezers

QB-6646  
Mouth Mirror

QB-6647

First Aid Examination Set of 3 Pieces



QB-6648

First Aid Examination Set of 3 Pieces



QB-6649

First Aid Examination Set of 3 Pieces



QB-6650

First Aid Examination Set of 3 Pieces





QB-6651  
Surgical  
11.5cm



QB-6652  
London College  
15cm



QB-6653  
Semken Taylor  
Str.



QB-6654  
Semken Taylor  
Cvd.



QB-6655  
Meriam  
16cm



QB-6656  
Dressing  
14cm



QB-6657  
Semken Taylor  
1x2 Straight



QB-6658  
Semken Taylor  
1x2 Curved



QB-6659  
Gilles  
15cm 1x2



QB-6660  
14cm 1x2



QB-6661  
Perry  
13cm



QB-6662  
Atraumatic  
15cm Cvd



QB-6663  
Atraumatic  
15cm Str.



QB-6664  
College with Lock  
15cm Cvd.



QB-6665  
Gerald  
17.5cm



QB-6666  
with lock  
15cm



QB-6667  
Atraumatic  
17.5cm

## TWEEZERS



QB-6668  
Adson-TC  
12cm



QB-6669  
Adson-TC  
12cm 1x2



QB-6670  
Adson-TC  
15cm 1x2



QB-6671  
Adson  
12cm



QB-6672  
Adson  
12cm 1x2



QB-6673  
Adson  
15cm 1x2



QB-6674  
Adson Brown  
12cm



QB-6675  
Crane Kaplan  
15cm



QB-6676  
Crane Kaplan  
15cm



QB-6677  
Crane Kaplan  
15cm



QB-6678  
Hemingway  
15cm



QB-6679  
Miller  
15cm Str.



QB-6680  
Miller  
15cm Cvd.



QB-6681  
Dressing  
Tweezer  
14.5cm



QB-6682  
College  
Curved  
Tip 15cm



QB-6683  
Maryam  
Tweezers  
15cm



QB-6684  
College  
Tweezer  
15cm




QB-6685  
Perry  
Tweezer  
12.5cm




QB-6686  
Self-locking  
Tweezer  
15cm

Single Ended Probe Solid Handle 4.5mm Length 150mm




Explorers-Single-End


Double-Ended Explorers



Probes Double Ended 6mm, length 170mm




Pointed Probes Double Ended Handle 6mm, length 170mm




EXPLORERS INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

PLUGGERS



## PERIODONTAL PROBES




## PERIODONTAL PROBES

### ENDODONTICS INSTRUMENTS

#### ANTERIOR





**GRACEY MINI FIVE**

Hollow Handle 9.5mm, Length 170mm Ceramic Coated Handle



QB-6777  
Fig. 3/4 Gracey 3M



QB-6778  
Fig. 5/6 Gracey 3M



QB-6779  
Fig. 7/8 Gracey 3M  
Buccal and lingual/  
Palatal surfaces



QB-6780  
Fig. 11/12 Gracey 3M  
Mesial



QB-6781  
Fig. 13/14 Gracey 3M  
Distal

QB-6776  
Fig. 1/2  
Gracey 3m



**GRACEY AFTER FIVE**

Hollow Handle 9.5mm, Length 170mm Ceramic Coated Handle



QB-6782  
Gracey 3/4 AF



QB-6783  
Gracey 5/6 AF



QB-6784  
Gracey 7/8 AF



QB-6785  
Gracey 11/12 AF



QB-6786  
Gracey 13/14 AF

QB-6787  
Gracey  
1/2 AF



**Gracey Curettes**




■ 8112 - 8141 - 8146  
 ■ 8150 - 8151  
 ■ 8152 - 8153  
 ■ 8154 - 8155



■ 8152 - 8141 - 8146  
 ■ 8150 - 8151  
 ■ 8152 - 8153  
 ■ 8154 - 8155

**Sickle Scalars**



QB-6797 Fig. U15/33 - Towner-Jaquette



QB-6798 Fig. U15/30 - Towner-Jaquette



QB-6799 Fig. H5/33



QB-6800 Fig. H6/H7



QB-6801 Fig. 30/33 - Jaquette



QB-6802 Fig. T2/T3 - Taylor



QB-6803 Fig. 204



QB-6804 Fig. 204S



QB-6805 Fig. 11/12 - MC - McCall  
Mesial



QB-6806 Fig. 13/14 - MC - McCall  
Distal



**Straight and curved sickle scalars**  
Feature: 2 cutting edges curved facial component tipped ends



**Universal curettes**  
are characterized by an angle of 90° between the terminal stem and the cutting component.



**Gracey curettes**  
are characterized by an angle of 70° between the terminal stem and the cutting component.

**Universal Curettes**



QB-6807 Fig. 13S/14S - MC - McCall



QB-6808 Fig. 17S/18S - MC - McCall  
Distal



QB-6809 Fig. M23



QB-6810 Fig. 204SD



QB-6811 Fig. YG7/8-Younger Good



QB-6812 Fig. L 1/2 - Langer



QB-6813 Fig. L 3/4 - Langer



QB-6814 Fig. L 5/6 - Langer



QB-6815 Fig. M 23A




QB-6816 Fig. M 23B



**Straight and curved sickle scalers**  
Feature: 2 cutting edges curved facial component tipped ends



**Universal curettes**  
are characterized by an angle of 90° between the terminal stem and the cutting component.



**Gracey curettes**  
are characterized by an angle of 70° between the terminal stem and the cutting component.

Periodontal knives





QB-6825  
17



QB-6826  
17D



QB-6827  
54A  
Mühlemann



QB-6828  
67



QB-6829  
37D



QB-6830  
8  
Wiland



QB-6831  
8  
Holz

0.9 MM



QB-6832  
9  
Holz

0.7 MM

Chisels

Diamond coated margin trimmers



QB-6833  
40P



QB-6834  
40A

Thanks to the diamond-coated hand instruments the margins are evenly levelled.

Margin trimmers



QB-6835  
1TH2 Hirt



QB-6836  
5TH6 Hirt



QB-6837  
3TH4 Hirt




QB-6838  
7TH8 Hirt



QB-6839  
62NS61



QB-6840  
37NS36



QB-6841  
66NS65



QB-6842  
79CI80 Black



QB-6843  
77CI78 Black



QB-6844  
279CI280 Black



QB-6845  
277CI278 Black



QB-6846  
179CI180 Black



QB-6847  
177CI178 Black

**CHISELS**



QB-6848  
1/2 Werdelstaedt



QB-6849  
3/4 Werdelstaedt



QB-6850  
7 Straight



QB-6851  
5/6 Straight



QB-6852  
7/10 Straight



QB-6853  
8/9 Binangle



QB-6854  
11/12 Binangle



QB-6855  
40/41 Binangle

**Tartar Remover Instruments Solid Handle Ø 6 mm**



QB-6856  
Fig. 1



QB-6857  
Fig. 1L  
Left



QB-6858  
Fig. 1R  
Right



QB-6859  
Fig. 2L  
left



QB-6860  
Fig. 2R  
Right



QB-6861



QB-6862



QB-6863

## BURNISHERS CONSERVATIVE TREATMENT



QB-6864  
WHITE 28



QB-6865  
WHITE 29




QB-6866  
WHITE 31



QB-6867  
WHITE 32




QB-6868  
WHITE 33



QB-6869  
WHITE 34



QB-6870  
WESTCOTT 21



QB-6871  
WESTCOTT 21B




QB-6872



QB-6873



QB-6874  
WHITE



QB-6875  
WHITE



QB-6876



QB-6877



QB-6878



QB-6879



QB-6880



QB-6881



QB-6882




QB-6883




QB-6884



QB-6885  
PLANATOR



QB-6886  
PLANATOR



## BURNISHERS CONSERVATIVE TREATMENT


### Ball Burnishers



### Hatches




## COMPOSITE INSTRUMENTS



## COMPOSITE INSTRUMENTS

### Spatulas





### Spatulas-Pluggers




## EXCAVATOR INSTRUMENTS


Water drop Our excavators have particularly sharp blades, free of burrs and perfectly rounded.








## COMPOSITE BALL PLUGGERS













QB-7005  
WILLIAMS



QB-7006  
WILLIAMS



QB-7007  
WILLIAMS



QB-7008



QB-7009



QB-7010



QB-7011



QB-7012



QB-7013




QB-7014

## COMPOSITE PLUGGERS INSTRUMENTS



## COMPOSITE PLUGGERS INSTRUMENTS



## Serrated Pluggers



|          |  |                  |          |
|----------|--|------------------|----------|
| Ø 0.9 MM |  | QB-7034<br>OP36M | Ø 1.2 MM |
| Ø 1.1 MM |  | QB-7035<br>OP36D | Ø 1.5 MM |
| Ø 1.6 MM |  | QB-7036<br>Op36  | Ø 2.2 MM |
| Ø 1.0 MM |  | QB-7037<br>Op52  | Ø 1.7 MM |
| Ø 2.2 MM |  | QB-7038<br>Zi51  | Ø 2.6 MM |
| Ø 1.8 MM |  | QB-7039<br>Op44  | Ø 2.2 MM |
| Ø 1.7 MM |  | QB-7040<br>1BR2  | Ø 2.6 MM |
| Ø 2.6 MM |  | QB-7041<br>2BR3  | Ø 3.6 MM |



## COMPOSITE MODELLING INSTRUMENTS



QB-7041 MF2C  
Fuzzi



QB-7043  
PJ4



QB-7044 MB3C  
Besek



QB-7045  
HZ3C



QB-7046  
CPRO3C



QB-7047 Mo1  
Monari



QB-7048 TH28C  
Hirt



QB-7049  
HZ25C



QB-7050  
Ck10



QB-7051  
CK28BC



QB-7052 Op37  
Carver



QB-7053  
28BR29

## Round excavators

### Cariology

- for the removal of softened, carious dentin
- for the removal of temporary fillings

### Surgery

- to sever the periodontal ligament prior to extraction
- for the removal of granulation tissue after extraction

Design: Round working end Curved cutting edge



Round Excavator, 1.0 mm  
QB-7054



Round Excavator, 1.5 mm  
QB-7055



Round Excavator, 2.0 mm  
QB-7056

## Triple-angled back action excavators



Triple-angled Back Action Excavator, 1.0 mm  
QB-7057



Triple-angled Back Action Excavator, 1.5 mm  
QB-7058



Triple-angled Back Action Excavator, 2.0 mm  
QB-7059



Triple-angled Back Action Excavator, 3.5 mm  
QB-7060




## FILLING INSTRUMENTS

COMPOSITE/PLASTIC FILLING INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

BLACK


PLUGGERS/CONDENSERS

Hollenback



MARQUETTE CONDENSER


MORTONSON



## FILLING INSTRUMENTS


RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM  
Back Action Filling Instruments

### Tanner Filling Instruments



QB-7074  
2T  
Tanner

QB-7075  
3T  
Tanner



QB-7076  
9/10  
Back Action

QB-7077  
11/12  
Back Action


QB-7078  
153  
English Pattern  
Serrated

QB-7079  
153  
English Pattern  
Serrated



QB-7073  
1T  
Tanner

### BALL



QB-7080  
18

QB-7081  
26/27S

QB-7082  
27/29

### BURNISHERS

### BEAVERTAIL



QB-7083  
2

QB-7084  
2  
Double-End

### ACORN



QB-7085  
21

QB-7086  
21B

### LADMORE



QB-7087  
2  
Ladmore

## CARVERS INSTRUMENTS

### P.K. THOMS WAXING INSTRUMENTS SOLID HANDLE 6mm




QB-7094  
Composite Instruments Blue  
Titanium with Cassettes



QB-7095  
P.k. Thomas Waxing Instruments Kit

### COMPOSITE INSTRUMENTS TITANIUM BLUE ENDODONTIC



## CARVERS INSTRUMENTS



QB-7102



QB-7103



QB-7104



QB-7105



QB-7106




Cassette For 5 Instruments

QB-7107

## FILLING INSTRUMENTS

FILLING INSTRUMENTS (Gold Pointed) SOLID HANDLE 6MM, LENGTH 170MM





QB-7118

QB-7119

QB-7120

QB-7121

QB-7122

QB-7123



QB-7124

QB-7125

QB-7126

## COMPOSITE FILLING INSTRUMENTS

### Composite Gold Tip Instruments



## MISCELLANEOUS PRODUCTS



QB-7138  
X Ray Film  
Hangers



QB-7139  
X Ray Film Hangers



QB-7140  
X Ray Film Hangers



QB-7141  
X Ray Film Hangers



QB-7142  
X Ray Film Hangers



QB-7143  
X Ray Film Hangers



QB-7144  
X Ray Film Hangers



QB-7145



QB-7146



QB-7147



Dissecting Set  
QB-7148 Set of 9 Pcs.  
QB-7149 Set of 11 Pcs.



QB-7150  
Napkin Holder with  
Plastic Chain

QB-7151  
Napkin holder with metal chain

### Key Rings




Mouth Props  
QB-7152  
Small

Mouth Props  
QB-7153  
Medium

Mouth Props  
QB-7154  
Large



QB-7155



QB-7156



QB-7157



QB-7158



QB-7159



QB-7160

## CEMENT SPATULAS



QB-7161  
00



QB-7162  
0



QB-7163  
1



QB-7164  
2



QB-7165  
3



QB-7166  
4



QB-7167  
3  
Weston



QB-7168  
1  
medium cement



QB-7169  
1



QB-7170  
Le Cron



QB-7171  
Le Cron



QB-7172  
Le Cron



QB-7173  
MSY



QB-7174  
MSY



QB-7175  
Le Cron



QB-7176  
Hylin



QB-7177  
Hylin



QB-7178  
Hylin



QB-7179  
Zahle



QB-7180  
Zahle



QB-7181  
Zahle



QB-7182  
Gritman

ALGinate, PLASTER SPATULAS, WAX KNIVES PLASTER FOR KNIVES & PLASTER SHEAR



QB-7183  
Wooden Handle  
Fig. 5  
21.5 cm



QB-7184  
Wooden Handle  
Fig. 7  
19.5 cm curved



QB-7185  
Wooden Handle  
Fig. 8  
20 cm curved



QB-7186  
Wooden Handle  
Fig. 9  
22.5 cm curved



QB-7187  
Hook Knives



QB-7188  
Fig. 1  
Lessmann  
12.5 cm



QB-7189  
Fig. 2  
Lessmann  
17 cm



QB-7190  
Fig. 1  
with small spoon  
12.5 cm



QB-7191  
Fig. 2  
with large spoon  
17 cm



QB-7192  
Bacolite Handle  
Wax Knives,  
Small



QB-7193  
Bacolite Handle  
Wax Knives,  
Large



QB-7194  
Wooden Handle  
Wax Knives,  
Small



QB-7195  
Wooden Handle  
Wax Knives,  
Large



QB-7196  
Plaster Knife



QB-7197  
Fig. 1  
17 cm



QB-7198  
Fig. 3  
16 cm



QB-7199  
Fig. 4  
16 cm



QB-7200  
Special Inox  
19.5 cm



**QB-7201**

**Technical Datas:**

Intercondylar distance  
condylar inclination  
Bennet angle  
Incisal table inclination  
Distance between mounting plates  
Weight  
Opening (item 5000)  
Opening (item 5010)

120 mm.  
fixed at 30°  
fixed at 7°  
10°  
85 mm.  
400 gr.  
125°  
180°



**QB-7202**

**Technical Datas:**

Intercondylar distance  
condylar inclination  
Bennet angle  
Incisal table inclination  
Distance between mounting plates  
Weight  
Opening

110 mm.  
fixed at 30°  
fixed at 7°  
flat  
85 mm.  
450 gr.  
180°



**QB-7203**

Stainless Steel Finishing

**QB-7204**

Royal Blue



**QB-7205**



**QB-7206**



**QB-7207**



**QB-7208**



**QB-7209**



**QB-7210**



**QB-7211**



**QB-7212**



**QB-7213**



**QB-7214**

**Nakawalni**  
Dental Laboratory  
Instrument

INSTRUMENTS  
FOR A HEALTHIER  
BRIGHTER SMILE

*Quality  
Instruments  
for a Brighter  
Tomorrow*



PROFESSIONALS  
TOOLS  
BRIGHTER  
SMILES  
WORLDWIDE

# TRUSTED DENTAL INSTRUMENT SOLUTIONS

MANUFACTURER | SUPPLIER | GLOBAL PARTNER

**Qbro Industries** is a trusted manufacturer and supplier of high-quality dental instruments, serving dental professionals worldwide. Our extensive range of diagnostic, periodontal, and restorative instruments is designed for precision, durability and clinical performance — helping you deliver better care and brighter smiles.

*Trusted  
today.  
Healthier  
tomorrow.*



**PRECISION  
CRAFTSMANSHIP**  
Engineered for  
Accuracy



**RELIABLE  
QUALITY**  
Consistent  
Performance



**CLINICAL  
INNOVATION**  
Supporting Better  
Outcomes



**WORLDWIDE  
SUPPORT**  
Your Partner  
Everywhere



MIRRORS



PROBES



SCALERS



CURETTES



TWEEZERS



SPATULAS



RESTORATIVE  
INSTRUMENTS

## Qbro Industries



Phone/WhatsApp:  
+4917620793624  
+923007705082  
+923039359066



Email: [info@qbroindustries.com](mailto:info@qbroindustries.com)



GLOBAL  
PARTNER  
IN DENTAL  
EXCELLENCE



QUALITY CARE  
WORLDWIDE